Gebruiksaanwijzing Operating Instructions Manual de instrucciones Mode d'emploi Istruzioni d'uso

G688053

CE0123

Rx only

medicon

Medicon eG · Gänssacker 15
78532 Tuttlingen · Germany
Tel.: +49 (0) 74 62 / 20 09-0
Fax: +49 (0) 74 62 / 20 09-50
E-Mail: sales@medicon.de
Internet: www.medicon.de

DE

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EN

Please read these instructions carefully before clinical application and keep them safe and close to hand. The instructions must be followed.

MEDICON implants are manufactured from pure titanium or Ti6Al4V alloy. Both materials are biocompatible, corrosionproof and nontoxic in the biological environment. They allow virtually artifact-free imaging through:

- conventional x-ray radiography
- computer tomography
- MRT (magnetic resonance tomography).

The surface is chemically passive and the material is antimagnetic. The implants are supplied in sterile condition and must be disinfected, cleaned and sterilized prior to use.

Table of contents:

- General information
- Purpose and proper use/indications
- Contraindications
- Possible side effects and complications
- Single use product
- CT and MRI instructions
- Implant selection
- Application / Handling of plates
- Application / Handling of screws
10. Instruments
11. Decontamination, cleaning and sterilization
12. Liability
13. Signs and Symbols

1. GENERAL INFORMATION
The implants are supplied sterile and must be disinfected, cleaned and sterilized prior to use. Please observe the following instructions. This will guarantee that the implant functions faultlessly and reliably.

2. PURPOSE AND PROPER USE/INDICATIONS
The MEDICON Ti-Epiplating CMS is intended for use in surgical interventions on the osseous structures of the skull for osseointegrated bone anchorage of craniofacial prosthesis and bone conduction hearing systems caused by craniofacial defects as:

- inherent malformation
- trauma
- after surgical tumor therapy
- reconstructive surgery

– in cases where other methods of treatment or devices were unsuccessful

The implants have been designed for osseointegration, not for replacing normal body structures or bearing the body weight in cases of insufficient osseointegration. The MEDICON Ti-Epiplating CMS may be used only by sufficiently trained and experienced surgeons.

3. CONTRAINDICATIONS

Patients that are unable to follow the instructions for postoperative care. This may be caused e.g. by:

- psychological/mental or neurological problems of the patient.
- Patients with inadequate bone tissue, or bone tissue of insufficient quality, and patients with circulatory disorders or latent infections.
- Extreme material sensitivity, i.e. foreign-body reaction of the patient. For this case, appropriate tests must be carried out under all circumstances prior to implantation (even when there is only a suspicion of such condition).
- Acute infections

4. POSSIBLE SIDE EFFECTS AND COMPLICATIONS

- Complaints, pain, abnormal sensations or palpability of the implant
- Extreme material sensitivity of the patient (foreign-body sensitivity) manifested by allergic reactions, etc.
- Heightened reaction of connective tissue in the region of the implant's position

– Inevitable in cases where other methods of treatment or devices were unsuccessful

– Insufficient bone formation, osteolysis, osteomyelitis, osteoporosis, inhibited revascularization or infections that can result in loosening, bending, rupture or breakage of the implants

– Delayed or insufficient osseointegration, which leads to a loosening of the implant

– Bone necrosis

– Inadequate healing. If the implants were not adapted properly, osseointegration can be delayed, is insufficient or does not occur at all

– Broken, bending, migration or loosening of the implant

– Reduced bone density due to stress shielding

Apart from the above and complications mentioned above, the surgical intervention can also lead to problems such as nerve injuries, infections, pain, etc., which are not directly attributable to the implant.

The possible complications are often the result of incorrect selection of the patient, lack of exercise or inadequate planning of the optimal implant position, rather than being caused by the implant itself. Extreme bending of the plates, or repeated bending at the same place, can lead to plate fractures. The implants can come loose if the screws are not tightened sufficiently.

5. SINGLE USE PRODUCT

Implants are designed and constructed for single use only on a single patient and may not be reused. An explanted implant must never be re-used. Even if the implant appears undamaged or functionally intact, indications of wear, small defects, and invisible stresses may exist. Since it is unclear what effect the forces and conditions within the body may have had on the stability, function, and material properties of an explanted implant, reimplantation incurs an unacceptable risk of early wearout or failure. The user is liable for non-compliance with the instructions for use.

The counter sink REF. 68.80.50 is designed and constructed for single use only on a single patient and may not be reused.

The counter sink must not be re-used, as it cannot be guaranteed that it remains sterile during preparation. Fine tissue remains that cannot be removed may collect in the fine structure of the counter sink. The user will be held liable if the instructions for Use are not observed.

6. CT AND MRI INSTRUCTIONS

The surgeon must inform the patient that he has to inform the responsible personnel about his implant when undergoing CT and MRI examinations. It is recommended to create an implant pass.

This can be used for information purposes when handed over to the responsible personnel by the patient. The implants pose no risk with regard to heating up and rotational behaviour during a CT or MRT with a field strength of up to 3 Tesla.

7. IMPLANT SELECTION

The surgeon is responsible for the correct selection of implants.

The following aspect is crucial for the correct selection of implants:

- Bone defect to be treated or intended bone repositioning
- Body weight of the patient

Incorrect implant selection can lead to premature implant loss and loosening, deformations or fractures of the implants. Among other factors, the success of the intervention depends on how the implants are handled. The plates must be bent carefully and gradually, avoiding extensive deformation. Damage and scratches reduce the stability of the product and result in early fatigue of the implant material. A product that has been used once must not be reused under any circumstances. Even if the product appears to be macroscopically undamaged, previous stresses on it may have led to quality loss and shortening of the lifespan of the product.

In cases of failed, delayed or insufficient osseointegration, the incidence of implant bending, implant failure or implant breakage cannot be excluded. A stress-free healing period of about 6 weeks is necessary and recommended. In cases of failed, delayed or insufficient osseointegration, the incidence of implant bending, implant failure or implant breakage cannot be excluded. A stress-free healing period of about 6 weeks is necessary and recommended. In cases of failed, delayed or insufficient osseointegration, the incidence of implant bending, implant failure or implant breakage cannot be excluded. A stress-free healing period of about 6 weeks is necessary and recommended.

We recommend using the "Emergency CMS screws" as standard screws for the Medicon Ti Epiplating CMS. There are no spare screws if a secure anchoring of the screws in the bone is not possible.

To ensure that the screwdriver blade is firmly seated in the screw head, axial force must be applied when inserting the screwdriver into the screw head. This ensures the correct axial alignment of screw and screwdriver, prevents slipping of, or damage to the screw head or the screwdriver blade. Application of excessive force in tightening the screws can lead to screw breakage during the intervention. Wear of the screwdriver blade

observation of this advice can put at risk any safe osseointegration. The patient must be instructed in turn, without delay, the operating surgeon of unusual changes at the operating site. The patient must be monitored carefully if there is a manifest change in the fixation region. The surgeon must be informed of any complications, e.g. implant failure, and discuss with the patient any necessary measures for further healing. The implants are not removed. An implant left in the body after complete healing. They are used for fixation of craniofacial prosthesis. It may act as a load-bearing element in active patients while thesurgion is obligated to control the correct fit of the implant and possibly make an exchange. For this reason, too, we accept liability for the implants only up to the period when the instruments must be inspected to make certain they are mechanically intact, not deformed, and fully functional. Do not use instruments that were found to be damaged, deformed or not fully functional. Exclusively instruments and accessories designated by Medicon are to be used in order to prevent any risks relating to the compatibility of the products.

11. DECONTAMINATION, CLEANING AND STERILIZATION
The following vJ/C/K prion-specific protective measure (z.B. for the instrument cleaning) is indicated.

In the event of diagnosis of a definite or probable case of vJ/C/K.
If it is not possible to use disposable products, the instrument used, which has been contaminated where contamination cannot be ruled out, must be disposed of as incinerated waste.

If prion contamination is suspected:
In the event of suspected prion contamination, incineration of the instrument is recommended according to the vJ/C/K task force final report.

If vJ/C/K is excluded:
Continued use after instrument processing is completed. Otherwise, the instrument, which has been contaminated or where contamination cannot be ruled out, must be disposed of as incinerated waste.

In the case of a non-identifiable vJ/C/K illness.
Even if nothing is known about the presence of a prion disease, two processing procedures should be used with at least partial efficacy against prions – e.g. mechanical alkaline cleaning combined with steam sterilization.

If mechanical alkaline cleaning or another cleaning procedure with proven efficacy against prions is not used and the medical devices in question are in contact with risk tissues (ZNS, eyes, lymphatic tissue), the RKI recommends a prolonged sterilization time of 18 minutes at 134°C.

Instruments made of stainless steel must not be placed in physiological saline solution (NaCl solution) as prolonged contact leads to corrosion such as pitting and stress-corrosion cracking. Only cleaned and disinfected implants may be sterilized.

Reprocessing restrictions
Frequent reprocessing has little effect on these instruments.

Location of use
Implants that have been contaminated by blood and/or secretions must not be placed back in the implant storage tray. They must be disposed of.

Preparation for decontamination
The instruments must be placed on appropriate instrument carriers so that they can be rinsed in the washing and disinfection device. The instrument carriers (e.g. sieve trays) must be such that subsequent cleaning in an ultrasonic device or in the washing and disinfection device is not impaired by any areas that are inaccessible to ultrasonics or rinsing.

Pre-cleaning
– Place implant storage tray for 15 minutes in ultrasonic bath at 40°C with 0.5% alkaline cleaner (Nedisher Mediconclean foa) and wash with ultra-sonic waves.

– Remove the instruments and rinse with cold, softened water.

A high contamination load in the ultrasonic cleaner requires the cleaning action and promotes the risk of corrosion. The cleaning solution must be renewed regularly according to the conditions of use. The criterion is visibly apparent soiling. In any case, a frequent change of bath is necessary, at least once a day. National guidelines must be observed.

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